

SALTER LABS®

Peak Flowmeters

Exciting Additions to Our Quality Products

Mini-Wright® Standard Range

Mini-Wright® AFS Low Range

PocketPeak® Universal Range

The health care provider can optimize accuracy and repeatability by selecting the model best suited to a patient's life style and PEF range.

All models are individually packaged and include:

- Re-usable mouthpiece
- Individual PEF daily diary
- Color keyed zone labels



Mini-Wright®
Standard Range
P/N 1235

The scalloped end piece design prevents a flowmeter from rolling off a counter top and onto the floor.

Large color keyed zone labels encourage a patient to develop a daily monitoring regimen. It also allows the care provider to customize the treatment regimen for each patient and circumstance.



Mini-Wright®
AFS Low Range
P/N 1236

Indicator arrows are recessed to prevent accidental breakage and are easily re-zeroed at the touch of a finger tip.

All Mini-Wright® PEF Meters labeled with blue or purple scales meet or exceed NAEP* published guidelines for accuracy, reproducibility and inter-device variability; as defined in the "statement for technical standards for Peak Flowmeters published in October 1992".

Ordering Information

Mini-Wright® Standard Range 60-880 L.P.M. 12/case P/N 1235

Mini-Wright® AFS Low Range 30-400 L.P.M. 12/case P/N 1236

Be sure to ask your Salter Labs® sales specialist about our full product line of medication nebulizers, treatment compressors, and MDI Chambers

Variety of uses. . . Home, Office, School, Clinic, Physician's Office



slc-62

SALTER LABS®

POCKETPEAK® for Asthma Management

Exciting Additions to Our Quality Products

When you compare the POCKETPEAK® with other peak flowmeters, compact size is not the only advantage. The POCKETPEAK® also excels in ease of use, performance and reliability.

POCKETPEAK® is so compact and easy to use it encourages patient compliance.



POCKETPEAK®
Universal Range
P/N 1230

POCKETPEAK's **ergonomic design** eliminates user confusion and misuse. The reusable mouthpiece attaches to the POCKETPEAK® for secure storage.

The ABS plastic housing and stainless steel vane are **dishwasher safe**. Cleaning the POCKETPEAK® and reusable mouthpiece is easy and worry free.

Universal Range allows the Health Care Professional to provide the most appropriate model for their patient.

POCKETPEAK® Flowmeter meets or exceeds all **NAEP* and NHLBI standards** and repeated testing over thousands of cycles demonstrates unsurpassed consistency, accuracy and repeatability of readings.

Ordering Information

POCKETPEAK® Universal Range	12/case	P/N 1230
Includes...Daily Record Diary		
Color Keyed Zone Labels		
Reusable Mouthpiece		

Specifications

Size: 3.4 inches x 3.2 inches
Weight: 1.4 Ounces
Range: 50-720 L.P.M.
Scale Increments: 10 L.P.M.
Repeatability: ± 2% L.P.M. at any reading

**The National Asthma Education Program (NAEP) guidelines require that a peak flowmeter be accurate over the full range ± 10% of the readings. Accuracy should be verified using an ATS Waveform generator. Reproducibility should be within 10 litres per minute or ± 5% of the readings, whichever is the larger. Inter-device variability should be ± 5%.*

Ref. National Asthma Education Program Statement on Technical Standards for Peak Flowmeters, January 1991

SALTER LABS®

100 W. Sycamore Road - Arvin, CA 93203
Phone: 661-854-3166
1-800-421-0024
1-800-235-4203
Fax: 661-854-3850
Toll Free Fax: 1-800-628-4690
Fax: 661-854-3850
E-mail: salterlabs@us.salterlabs.com
www.salterlabs.com

CE
0482

ISO 13485