


SALTER LABS®

Pocket Chamber® A Compact Valved Chamber for MDI

Exciting Addition to Our Quality Products

Unique, compact, patented holding chamber device for enhancing medication delivery to patients with treatment regimens utilizing metered dose inhalers (MDI)



Pocket Chamber®

A compact
valved chamber
for MDI
P/N 1240

Features

- Uni-body construction
- Produced from Anti-static material
- Rounded shoulder design
- Patented one way valve (chamber)
- Unique grid plate (spacer)
- Whistle alarm
- Unique design
- Latex free
- Sensitive one-way valve

Benefits

- Durability, easy to clean
- Clear, low electrostatic charge provides consistent dose delivery
- Re-shapes aerosol plume, improves particle distribution
- Very low resistance allows maximum performance
- Decelerates plume, improves particle distribution in small body
- Alerts patient when optimum inspiratory rate is exceeded
- Pocket-sized for patient convenience
- Eliminates potential patient allergies to latex
- Prevents inadvertent exhalation into the chamber


Variety of uses. . . Home, Office, School, Clinic, Physician's Office



slc-30

Pocket Chamber®

- Delivered Dose Equivalency at a Competitive Cost
- Delivered Dose Equivalency in a More Convenient Size



- Compact 110 cc size, pocketable easy to carry and use
- Tethered mouthpiece lintcap
- A universal style PVC end piece conforms to all popular MDI actuators
- Maximize mg/ml. per available dose
- Instructions on body, improves patient compliance

Ordering Information	units per case	part #
Pocket Chamber (one way valve design)	10 or 50	1240



100 W. Sycamore Road - Arvin, CA 93203
Phone: 661-854-3166
1-800-421-0024
1-800-235-4203
Fax: 661-854-3850
Toll Free Fax: 1-800-628-4690
Fax: 661-854-3850
E-mail: salterlabs@us.salterlabs.com
www.salterlabs.com