


SALTER LABS®


Oxygen Delivery Products

Exciting Additions to Our Quality Products




CGA870 Oxygen Regulator
P/N #1275

Yoke type CGA 870 Regulator. Contents Gauge with eleven preset dial settings for delivery flow rates of 1/8, 1/4, 1/2, 1.5, 2, 2.5, 3, 4, 5, 6, or 8 LPM, rugged, reliable, cost effective. Long life, trouble free



CGA870 Oxygen Regulator
P/N #1276

All brass chrome plated, high performance 0-8 LPM yoke type CGA870 connection. Contents and LPM gauges. Top of the line regulator.



CGA540 Oxygen Regulator
P/N #1278

All brass chrome plated, 0-8 LPM contents and LPM gauge CGA540 hex nut connection. A top of the line quality and performance regulator.



CGA540 Oxygen Regulator
P/N #1277

CGA 540 H T connection, contents gauge. Eleven pre-set dial settings delivery flow rates of 1/8, 1/4, 1/2, 1, 1.5, 2, 3, 4, 5, 6, or 8 LPM. Easy to operate, rugged, cost effective regulator, Long life, trouble free.



Litemeters
P/N #1271,
#1272,
#1273

Convenient pocket size, designed to validate flow rates for compressors, regulators, flowmeters. Three flow rate ranges 0-2.5 LPM, 0-8 LPM and 6-15 LPM. An inexpensive way to assure your equipment is operating properly.

Variety of uses. . . Home, Hospital, Clinic, Physician's Office, Skilled Nursing Facility, Patient Transport





D-E Style Cylinder Cart
P/N #1265

**D-E Style Cylinder
Cart**
P/N #1265

Lightweight D or E cylinder cart. Constructed from Heavy wall 7/8" and 3/4" Aluminum tube stock. Mounted on 6" rubber wheels with bird cage style cylinder holder and nylon safety strap. A spring loaded (3) position adjustable handle makes this a convenient travel cylinder cart.



Cylinder Wrench
P/N #1246

Cylinder Wrench


Packaged 10.
P/N# 1246
Packaged 10 with chain.
P/N# 1246C

Salter Cylinder Wrench is a versatile tool for opening and closing post-valve cylinders. In addition, the comfortable shape aids in connecting "tee" bar style CGA870 regulators to the cylinder. It also has a slot for storage of large and small sealing washers.

Inlet Seal Washer
(Viton rubber/brass)
Packaged 5
P/N# 1248

Open End Wrench
P/N #1247

Easy to use open end wrench. Fits hex nut on regulators used with CGA540 oxygen cylinders. Constructed of durable metal.
P/N #1247 (not pictured)



Cylinder Wrench
P/N #1249C

Cylinder Wrench
Packaged 50.
P/N# 1249C

Basic Cylinder Wrench with security cord
CGA870



100 W. Sycamore Road - Arvin, CA 93203
Phone: 661-854-3166
1-800-421-0024
1-800-235-4203
Fax: 661-854-3850
Toll Free Fax: 1-800-628-4690
Fax: 661-854-3850
E-mail: salterlabs@us.salterlabs.com
www.salterlabs.com

CE
0482

ISO 13485